

# Accessory Compatibility Guide

|                                                                                                             |                                                                                   |                                                                                    |                                                                                     |                                                                                     |
|-------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| <div>VITALITY<sup>AI</sup></div> <div>AUDIBEL   Signature<br/>SERIES</div> <div>INTRIGUE<sup>AI</sup></div> | StarLink<br>Charger 2.0                                                           | StarLink Premium<br>Charger 2.0                                                    | StarLink Custom<br>Charger 2.0                                                      | StarLink Premium<br>Custom Charger 2.0                                              |
|                                                                                                             |  |  |  |  |
|                                                                                                             | RIC RT: P00001445<br>mRIC R: P00009119                                            | RIC RT: P00009123<br>mRIC R: P00009122                                             | ITC R, ITE R, CIC R NW:<br>P00001455                                                | ITC R, ITE R, CIC R NW:<br>P00003628                                                |
|                                                                                                             | Default*                                                                          |                                                                                    |                                                                                     |                                                                                     |
|                                                                                                             | Onboard Battery                                                                   |                                                                                    |                                                                                     |                                                                                     |
|                                                                                                             | Turbo Charge                                                                      |                                                                                    |                                                                                     |                                                                                     |

|                             |                                                                                    |                                                                                     |
|-----------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| <div>Arc<sup>AI</sup></div> | Standard and<br>Custom Charger                                                     | Standard and Custom Premium<br>Charger                                              |
|                             |  |  |
|                             | RIC R, BTE R: 43084-100<br>ITC R, ITE R: 43085-100                                 | RIC R, BTE R: 43014-100<br>ITC R, ITE R: 43041-100                                  |
|                             | Default*                                                                           |                                                                                     |
|                             | Onboard Battery                                                                    |                                                                                     |
|                             | Turbo Charge                                                                       |                                                                                     |

# Accessory Compatibility Guide

|                      |                        |                                                                                   |                                                                                   |                                                                                   |                                                                                     |                                                                                     |                                                                                     |
|----------------------|------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| WIRELESS ACCESSORIES |                        | StarLink Edge<br>TV Streamer                                                      | StarLink Remote<br>Control 2.0                                                    | StarLink Mini<br>Remote Microphone                                                | StarLink Remote<br>Microphone +                                                     | StarLink Table<br>Microphone                                                        | StarLink<br>TV Streamer                                                             |
|                      |                        |  |  |  |  |  |  |
|                      |                        | P00004316                                                                         | P00002347                                                                         | 43121-000                                                                         | 43120-000                                                                           | 43123-000                                                                           | 43119-001                                                                           |
|                      | VITALITY <sup>AI</sup> | ●                                                                                 | ●                                                                                 |                                                                                   |                                                                                     |                                                                                     |                                                                                     |
|                      | INTRIGUE <sup>AI</sup> |                                                                                   | ●                                                                                 | ●                                                                                 | ●                                                                                   | ●                                                                                   | ●                                                                                   |
|                      | Arc <sup>AI</sup>      |                                                                                   | ●                                                                                 | ●                                                                                 | ●                                                                                   | ●                                                                                   | ●                                                                                   |

|                        |                                                                                       |                                                                                       |                                                                                       |
|------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|                        | APPS                                                                                  |                                                                                       |                                                                                       |
|                        | My<br>Audibel                                                                         | Thrive                                                                                | Hear<br>Share                                                                         |
|                        |  |  |  |
| VITALITY <sup>AI</sup> | ●                                                                                     |                                                                                       | ●                                                                                     |
| INTRIGUE <sup>AI</sup> | ●                                                                                     |                                                                                       | ●                                                                                     |
| Arc <sup>AI</sup>      |                                                                                       | ●                                                                                     | ●                                                                                     |